

GOAT ANTI-SYNAPTOTAGMIN-10 ANTIBODY

SKU: EB11611



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	sytx SYT10 synaptotagmin-10 synaptotagmin X synaptotagmin 10 OTTHUMP00000240621 MGC119437 MGC119436
Accession ID	NP_945343.1
Blocking Peptide	EBP11611
Immunogen	Peptide with sequence CRTGLDAEGLGRDH, from the internal region (near C Terminus) of the protein sequence according to NP_945343.1.
Peptide Sequence	CRTGLDAEGLGRDH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig
Human Gene ID	341359
Mouse Gene ID	54526
Rat Gene ID	60567
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)