


GOAT ANTI-SYNAPTOTAGMIN-10 (MOUSE) ANTIBODY

SKU: EB11610



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	sytx Syt10 synaptotagmin-10 synaptotagmin X synaptotagmin 10
Usage	Flow Cytometry: Flow cytometric analysis of NIH3T3 cells. Recommended concentration:
Summary	10ug/ml.
Accession ID	NP_061273.1
Blocking Peptide	EBP11610
Immunogen	Peptide with sequence RKEDGVSSLCQK, from the N Terminus of the protein sequence according to NP_061273.1.
Peptide Sequence	RKEDGVSSLCQK

Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat, Dog, Pig, Cow
Reactive Species	Mouse
Mouse Gene ID	54526
Rat Gene ID	60567
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA, FC

GALLERY IMAGES

