

GOAT ANTI-SYNTAXIN 11 ANTIBODY

SKU: EB12336

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	syntaxin-11 HPLH4 HLH4 FHL4 syntaxin 11 STX11
Accession ID	NP_003755.2
Blocking Peptide	EBP12336
Immunogen	Peptide with sequence C-KQADTLNVIELNVQK, from the internal region of the protein sequence according to NP_003755.2.
Peptide Sequence	C-KQADTLNVIELNVQK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Dog
Reactive Species	Human
Human Gene ID	8676
Rat Gene ID	292483
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 36kDa band observed in Human Placenta lysates (calculated MW of 33.2kDa according to NP_003755.2). This molecular weight is routinely observed by other sources. Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

