

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-SYNTAXIN 11 ANTIBODY**

**SKU:** EB12336



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** syntaxin-11|HPLH4|HLH4|FHL4|syntaxin 11|STX11

Names

**Accession** 

NP\_003755.2

**Blocking** 

EBP12336 **Peptide** 

Peptide with sequence C-KQADTLNVIELNVQK, from the internal region of the protein sequence according to **Immunogen** 

NP 003755.2.

**Peptide** Sequence

C-KQADTLNVIELNVQK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

**Species** 

Human, Rat, Dog

Reactive

**Species** 

Human

Human

8676 Gene ID

**Rat Gene ID** 292483

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

Approx 36kDa band observed in Human Placenta lysates (calculated MW of 33.2kDa according to Western

NP\_003755.2). This molecular weight is routinely observed by other sources. Recommended concentration: Blot

 $1-3\mu g/ml$ .

**Application** 

Pep-ELISA, WB **Type** 

## **GALLERY IMAGES**



 $\textbf{Email:} \ \underline{customerservice@vectorlabs.com}$ 

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa