

Telephone: (650) 697-3600

GOAT ANTI-SYNTROPHIN GAMMA 1 ANTIBODY

SKU: EB06261

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias SNTG1|SYN4|G1SYN|syntrophin, gamma 1|gamma1-syntrophin|syntrophin 4 gamma1-syntrophin|syntrophin 4

Names

Accession NP 061840.1

Blocking

EBP06261 Peptide

Peptide with sequence C-ASSATTSKAKYTT, from the C Terminus of the protein sequence according to Immunogen

NP 061840.1.

Peptide

C-ASSATTSKAKYTT Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

54212 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:8000.

Limit

Preliminary experiments gave no signal but low background in Human Testis and Brain extracts at up to 1µg/ml. Western

We would appreciate any feedback from people in the field - have any results been reported with other Blot

antibodies/lysates?

Application

Pep-ELISA **Type**