

GOAT ANTI-SYNTROPHIN GAMMA 1 ANTIBODY

SKU: EB06261



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SNTG1 SYN4 G1SYN syntrophin, gamma 1 gamma1-syntrophin syntrophin 4 gamma1-syntrophin syntrophin 4
Accession ID	NP_061840.1
Blocking Peptide	EBP06261
Immunogen	Peptide with sequence C-ASSATTSKAKYTT, from the C Terminus of the protein sequence according to NP_061840.1.
Peptide Sequence	C-ASSATTSKAKYTT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	54212
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Preliminary experiments gave no signal but low background in Human Testis and Brain extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)