

GOAT ANTI-SYT7 ANTIBODY

SKU: EB10332

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Syt7 synaptotagmin VII synaptotagmin 7 B230112P13Rik A1851541
Accession ID	NP_004191.2
Blocking Peptide	EBP10332
Immunogen	Peptide with sequence C-DIPTEKLRET, from the internal region of the protein sequence according to NP_004191.2.
Peptide Sequence	C-DIPTEKLRET
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	9066
Mouse Gene ID	54525
Rat Gene ID	59267
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)