



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-SYT7 ANTIBODY**

**SKU:** EB10332

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Syt7|synaptotagmin VII|synaptotagmin 7|B230112P13Rik|AI851541

**Names** 

Accession NP 004191.2 ID

**Blocking** 

EBP10332 **Peptide** 

Immunogen Peptide with sequence C-DIPTEKLRET, from the internal region of the protein sequence according to NP 004191.2.

**Peptide** C-DIPTEKLRET Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** Refrigerated

Instructions

**Predicted** Human, Mouse, Rat, Dog, Cow

Human

9066 Gene ID

Mouse 54525 **Gene ID** 

Rat Gene ID 59267

**Product** Grade

**Species** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

**Application** Pep-ELISA **Type**