

GOAT ANTI-TAF1 ANTIBODY

SKU: EB05600



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	TATA box binding protein (TBP)-associated factor, RNA polymerase II, A, 250kD cell cycle, G1 phase defect, (transcription factor TFIID p250 polypeptide) complementation of cell cycle block, G1-to-S TAFII250 TAF2A NSCL2 P250 CCGS CCG1 BA2R TAF1 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 250 kD TAF1
Accession ID	NP_004597; NP_620278
Blocking Peptide	EBP05600
Immunogen	Peptide with sequence GPGCDLLLRTAAT, from the N Terminus of the protein sequence according to NP_004597; NP_620278.
Product Comments	This antibody is expected to recognise both human isoforms.
Peptide Sequence	GPGCDLLLRTAAT
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	6872
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 250kDa band seen in HeLa lysate (predicted MW of 208kDa according to NP_004597). Recommended for use at 0.5-2µg/ml.