



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-TAF4 ANTIBODY**

**SKU:** EB05596

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

TATA box binding protein (TBP)-associated factor, 135 kD|TAF4A RNA polymerase II|TATA box binding protein

(TBP)-associated factor, RNA polymerase II, C1, 130kD|transcription initiation factor TFIID 135 kD

Alias subunit|TAFII135|TAFII130|TAF2C1|TAF4A|TAF2C|TAF4: TAF4 RNA polymerase II, TATA box binding protein **Names** 

(TBP)-associated factor, 135 kD

**Accession** 

NP\_003176

**Blocking** 

EBP05596 **Peptide** 

Peptide with sequence C-TSHSLLLYKAFLK, from the C Terminus of the protein sequence according to **Immunogen** 

NP\_003176.

**Peptide** 

C-TSHSLLLYKAFLK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human, Mouse **Species** 

Reactive

Human **Species** 

Human Gene ID

6874

Mouse

117169

**Gene ID Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

Western Western blot: Approx 125kDa band seen in Human Brain lysate (Predicted MW of 119kDa according to

Blot NP 003176). Recommended for use at 1-3µg/ml.