



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-TAIP-12 / FAM130A1 ANTIBODY

SKU: EB07625



250kDa
150kDa
100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	family with sequence similarity 130, member A1 BC035295 cDNA sequence BC035295 MGC28165 TAIP-12 TGF-beta induced apoptosis protein 12 C12orf2 C12orf22 FLJ25576
Accession ID	NP_700456.1
Blocking Peptide	EBP07625
Immunogen	Peptide with sequence C-KRQKQLRRKNVRFDQ, from the internal region of the protein sequence according to NP_700456.1.
Peptide Sequence	C-KRQKQLRRKNVRFDQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human, Mouse
Human Gene ID	81566
Mouse Gene ID	207785
Rat Gene ID	315308
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 60kDa band observed in Human Heart and Mouse Heart lysates (calculated MW of 58.5kDa according to NP_700456.1). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB



DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES





www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)