

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TAIP-12 / FAM130A1 ANTIBODY

SKU: EB07625







250kDa 150kDa 100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

family with sequence similarity 130, member A1|BC035295|cDNA sequence

Alias Names

BC035295|MGC28165|TAIP-12|TGF-beta induced apoptosis protein 12|C12orf2|C12orf22|FLJ25576

Accession

ID

NP 700456.1

Blocking

Peptide

EBP07625

Immunogen

Peptide with sequence C-KRQKQLRRKNVRFDQ, from the internal region of the protein sequence according to

NP 700456.1.

Peptide

Sequence

C-KRQKQLRRKNVRFDQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Rat, Dog

Reactive

Species

Human, Mouse

207785

Human

81566 Gene ID

Mouse

Gene ID

Rat Gene ID 315308

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection Limit

Antibody detection limit dilution 1:64000.

Western

Blot

Approx 60kDa band observed in Human Heart and Mouse Heart lysates (calculated MW of 58.5kDa according

to NP_700456.1). Recommended concentration: 0.1-0.3µg/ml.

Application

Type

Pep-ELISA, WB



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa