

GOAT ANTI-TAPBPL ANTIBODY

SKU: EB05638



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	TAP binding protein-like TAPASIN-RELATED PROTEIN hypothetical protein FLJ10143 FLJ10143 TAPBPR TAP binding protein related TAPBP-R TAPBPL
Accession ID	NP_060479.3
Blocking Peptide	EBP05638
Immunogen	Peptide with sequence C-EDRTARVSQPS, from the C Terminus of the protein sequence according to NP_060479.3.
Peptide Sequence	C-EDRTARVSQPS
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	55080
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments gave an approx 60-65kDa band in Human Ovary and Placenta lysates at 1ug/ml. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 54kDa according to NP_060479. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

