



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TAPBPL ANTIBODY

SKU: EB05638

SPECIFICATIONS

Unit Size

100 µg

Synonym / Alias

TAP binding protein-like|TAPASIN-RELATED PROTEIN|hypothetical protein FLJ10143|FLJ10143|TAPBPR|TAP binding

protein related|TAPBP-R|TAPBPL Names

Accession

ID

NP_060479.3

Blocking

EBP05638 **Peptide**

Immunogen Peptide with sequence C-EDRTARVSQPS, from the C Terminus of the protein sequence according to NP_060479.3.

Peptide

Sequence

C-EDRTARVSQPS

Shipping

Instructions Refrigerated

Human

55080

Predicted

Species

Human

Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection

Western

Antibody detection limit dilution 1:64000.

Limit

Blot

Preliminary experiments gave an approx 60-65kDa band in Human Ovary and Placenta lysates at 1ug/ml. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 54kDda according to NP_060479. We would appreciate any feedback from people in the field have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms

been reported?

Application

Type

Pep-ELISA