Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TAS1R3 ANTIBODY

SKU: EB08728

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias T1R3|taste receptor, type 1, member 3|TAS1R3

Names

Accession NP 689414.1

Blocking

EBP08728 **Peptide**

Peptide with sequence C-GQNDGNTGNQGKHE, from the internal region of the protein sequence according to Immunogen

NP 689414.1.

Peptide C-GQNDGNTGNQGKHE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Human

Shipping Refrigerated Instructions

Predicted

Species

Human 83756 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this Blot

product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.

Application

Pep-ELISA **Type**