

GOAT ANTI-TAS1R3 ANTIBODY

SKU: EB08728



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	T1R3 taste receptor, type 1, member 3 TAS1R3
Accession ID	NP_689414.1
Blocking Peptide	EBP08728
Immunogen	Peptide with sequence C-GQNDGNTGNQGKHE, from the internal region of the protein sequence according to NP_689414.1.
Peptide Sequence	C-GQNDGNTGNQGKHE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	83756
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)