



# GOAT ANTI-TAS1R3 / T1R3 ANTIBODY

**SKU:** EB08223



---

## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias</b>	T1R3 taste receptor, type 1, member 3 TAS1R3
<b>Names</b>	
<b>Accession ID</b>	NP_689414.1
<b>Blocking Peptide</b>	EBP08223
<b>Immunogen</b>	Peptide with sequence C-QNDGNTGNQGKHE, from the C Terminus of the protein sequence according to NP_689414.1.
<b>Peptide Sequence</b>	C-QNDGNTGNQGKHE
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Human Gene ID</b>	83756
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:2000.
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS



[www.everestbiotech.com](http://www.everestbiotech.com)

Email: [customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

- [Data Sheet](#)