



Telephone: (650) 697-3600

GOAT ANTI-TAS1R3 / T1R3 ANTIBODY

SKU: EB08223

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias T1R3|taste receptor, type 1, member 3|TAS1R3

Names

Accession NP 689414.1 ID

Blocking

EBP08223 Peptide

Immunogen | NP_689414.1. Peptide with sequence C-QNDGNTGNQGKHE, from the C Terminus of the protein sequence according to

Peptide C-QNDGNTGNQGKHE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated Instructions

Predicted Human

Species Human

83756 **Gene ID**

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA Detection

Antibody detection limit dilution 1:2000. Limit

Application Pep-ELISA **Type**