

GOAT ANTI-TBC1D4 ANTIBODY

SKU: EB06099



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Acrg embryonic lethality minimal region ortholog DKFZp779C0666 AS160 RP11-159J2.3 AS160 Acrg embryonic lethality (mouse) minimal region ortholog TBC (Tre-2, BUB2, CDC16) domain-containing protein KIAA0603 gene product KIAA0603 TBC1 domain family, member 4 TBC1D4
Accession ID	NP_055647.2
Blocking Peptide	EBP06099
Immunogen	Peptide with sequence CNPNKAKIGNKP, from the C Terminus of the protein sequence according to NP_055647.2.
Peptide Sequence	CNPNNKAKIGNKP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	9882
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Human Pancreas and Kidney extracts at up to 1.5µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?