



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TBC1D4 ANTIBODY

SKU: EB06099

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Acrg embryonic lethality minimal region ortholog|DKFZp779C0666|AS160|RP11-159|2.3|AS160|Acrg embryonic

lethality (mouse) minimal region ortholog|TBC (Tre-2, BUB2, CDC16) domain-containing protein|KIAA0603 gene Alias

Names product|KIAA0603|TBC1 domain family, member 4|TBC1D4

Accession

NP 055647.2 ID

Blocking

EBP06099 **Peptide**

Immunogen Peptide with sequence CNPNNKAKIGNKP, from the C Terminus of the protein sequence according to NP 055647.2.

Peptide

CNPNNKAKIGNKP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human

Species

Human

9882 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection

Western

Antibody detection limit dilution 1:16000.

Limit

Blot

Western Blot: Preliminary experiments gave no signal but low background in Human Pancreas and Kidney extracts at up to 1.5µg/ml. We would appreciate any feedback from people in the field - have any results been reported

with other antibodies/lysates?