

GOAT ANTI-TBC1D9 (AA444-458) ANTIBODY

SKU: EB12302



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TBC1 domain family member 9A TBC1 domain family member 9 MDR1 TBC1 domain family, member 9 (with GRAM domain) TBC1D9
Accession ID	NP_055945.2
Blocking Peptide	EBP12302
Immunogen	Peptide with sequence C-DGERQFNLNGNSVPT, from the internal region of the protein sequence according to NP_055945.2.
Peptide Sequence	C-DGERQFNLNGNSVPT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Pig, Cow
Human Gene ID	23158
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)