Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TCEAL1 ANTIBODY

SKU: EB09260

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

OTTHUMP00000024295|OTTHUMP00000024294|OTTHUMP00000024293|pp21|p21|SIIR|transcription elongation **Alias**

factor A (SII)-like 1|TCEAL1 **Names**

Accession

NP 004771.2

Blocking

EBP09260 **Peptide**

Peptide with sequence DKPRKENEEEPQ-C, from the N Terminus of the protein sequence according to Immunogen

NP 004771.2.

Product

Reported variants represent identical protein: NP_001006641.1, NP_001006640.1, NP_004771.2. Comments

Peptide

DKPRKENEEEPQ-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Species

Human

Human

9338 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Application Pep-ELISA **Type**

> Goat Anti-TCEAL1 Antibody https://everestbiotech.com/products/goat-anti-tceal1-antibody/