

GOAT ANTI-TCEAL1 ANTIBODY

SKU: EB09260



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	OTTHUMP00000024295 OTTHUMP00000024294 OTTHUMP00000024293 pp21 p21 SIIR transcription elongation factor A (SII)-like 1 TCEAL1
Accession ID	NP_004771.2
Blocking Peptide	EBP09260
Immunogen	Peptide with sequence DKPRKENEEEPQ-C, from the N Terminus of the protein sequence according to NP_004771.2.
Product Comments	Reported variants represent identical protein: NP_001006641.1, NP_001006640.1, NP_004771.2.
Peptide Sequence	DKPRKENEEEPQ-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	9338
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)