

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TCF4 / ITF2 ANTIBODY

SKU: EB05283

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

immunoglobulin transcription factor 2|MGC149724|MGC149723|bHLHb19|Transcription factor-4

Alias (immunoglobulin transcription factor-2)|SEF2-1B|SEF2-1A|SEF2-1|SEF2|ITF2|E2-2|transcription factor 4|TCF4 **Names**

Accession

NP_001077431.1; NP_003190.1

Peptide with sequence C-HPGMGDASNHMGQM, from the C Terminus of the protein sequence according to **Immunogen**

NP 001077431.1; NP 003190.1.

Product

This antibody is expected to recognize both reported isoforms (NP_001077431.1and NP_003190.1). Comments

Peptide

Sequence

C-HPGMGDASNHMGQM

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

using the immunizing peptide. Method

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Dog

Reactive

Human **Species**

Human

6925 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Type

Western Blot: Approx 60kDa band observed in Human Placenta lysates (calculated MW of 71.3kDa according Western

Blot to NP 003190.1). Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600