

GOAT ANTI-TCHHL1 ANTIBODY

SKU: EB08714



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	intermediate filament-associated protein S100 calcium binding protein A17 basalin TCHHL1 S100A17 trichohyalin-like 1 TCHHL1
Accession ID	NP_001008536.1
Blocking Peptide	EBP08714
Immunogen	Peptide with sequence C-QALEDKQGHPQRE, from the internal region of the protein sequence according to NP_001008536.1.
Peptide Sequence	C-QALEDKQGHPQRE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	126637
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Approx. 90kDa band observed in Human Skin lysates (calculated MW of 99.3kDa according to NP_001008536.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

