



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TCHHL1 ANTIBODY

SKU: EB08714



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

intermediate filament-associated protein|S100 calcium binding protein Alias

A17|basalin|THHL1|S100A17|trichohyalin-like 1|TCHHL1 Names

Accession

NP_001008536.1 ID

Blocking

EBP08714 **Peptide**

Peptide with sequence C-QALEDKQGHPQRE, from the internal region of the protein sequence according to **Immunogen**

NP_001008536.1.

Peptide

C-QALEDKQGHPQRE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

126637 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:4000.

Limit

Western Approx. 90kDa band observed in Human Skin lysates (calculated MW of 99.3kDa according to

Blot NP_001008536.1). Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa