

GOAT ANTI-TCL1 (MOUSE) ANTIBODY

SKU: EB07721



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	T-cell lymphoma breakpoint 1 Tcl1
Accession ID	NP_033363.1; NP_001276397.1
Blocking Peptide	EBP07721
Immunogen	Peptide with sequence QLYPKDRYRSCD, from the internal region of the protein sequence according to NP_033363.1; NP_001276397.1.
Peptide Sequence	QLYPKDRYRSCD
Shipping Instructions	Refrigerated
Predicted Species	Mouse
Mouse Gene ID	21432
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments gave an approx 28kDa band in Rat Ovary lysates after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 14.1kDa according to Mouse NP_033363.1. The 28kDa band was successfully blocked by incubation with the immunizing peptide.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

