

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-TCL1 (MOUSE) ANTIBODY**

**SKU:** EB07721

×

## **SPECIFICATIONS**

Unit Size 100 μg

Synonym /

Alias T-cell lymphoma breakpoint 1|Tcl1

Names

ID

Accession

NP\_033363.1; NP\_001276397.1

Blocking EBP07721 Peptide

Peptide with sequence QLYPKDRYRSCD, from the internal region of the protein sequence according to

NP\_033363.1; NP\_001276397.1.

Peptide QLYPKDRYRSCD Sequence

Shipping

Instructions Refrigerated Predicted ...

Species Mouse
Mouse

Gene ID

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Preliminary experiments gave an approx 28kDa band in Rat Ovary lysates after 0.1µg/ml antibody staining. Please

Western

Blot

Preliminary experiments gave an approx 28kDa band in Rat Ovary lysates after 0.1µg/ml antibody staining. Please
note that currently we cannot find an explanation in the literature for the band we observe given the calculated
size of 14.1kDa according to Mouse NP\_033363.1. The 28kDa band was successfully blocked by incubation with