

GOAT ANTI-TCL1A ANTIBODY

SKU: EB07720



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TCL1A T-cell leukemia/lymphoma 1A TCL1 T-cell lymphoma-1 T-cell lymphoma-1A
Accession ID	NP_068801.1
Blocking Peptide	EBP07720
Immunogen	Peptide with sequence C-QLYPDGRYRSSD, from the internal region of the protein sequence according to NP_068801.1.
Peptide Sequence	C-QLYPDGRYRSSD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8115
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 16kDa band observed in lysates of cell lines Daudi and Jurkat (calculated MW of 13.5kDa according to NP_068801.1). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

