

GOAT ANTI-TEB4 / MARCH-VI (AA 701 TO 712) ANTIBODY

SKU: EB08294



250kDa
150kDa
100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	membrane-associated RING-CH protein VI RNF176 KIAA0597 membrane-associated ring finger (C3HC4) 6 MARCH6 TEB4 MARCH-VI
Accession ID	NP_005876.2
Blocking Peptide	EBP08294
Immunogen	Peptide with sequence C-QGRRVIFQKVKE, from the internal region of the protein sequence according to NP_005876.2.
Peptide Sequence	C-QGRRVIFQKVKE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	10299
Mouse Gene ID	223455
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 100kDa band observed in Human Lung lysates (calculated MW of 103kDa according to NP_005876.2). Recommended concentration: 0.2-0.5µg/ml. An additional band of unknown identity was also consistently observed at 18kDa. This band was successfully blocked by incubation with the immunising peptide.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



