


GOAT ANTI-TENSIN 1 (AA902-914) ANTIBODY

SKU: EB11098



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias DKFZp586K0617| matrix-remodelling associated 6| Matrix-remodelling-associated protein 6| MGC88584| MST091| MST122| MST127| MSTP122| MSTP127| MXRA6| OTTHUMP00000207627| OTTHUMP00000207628|

Names tensin 1| tensin-1| TNS|TNS1

Accession ID NP_072174.3

Blocking Peptide EBP11098

Immunogen Peptide with sequence C-QYENQSPEATSPR, from the internal region of the protein sequence according to NP_072174.3.

Peptide Sequence C-QYENQSPEATSPR

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat, Cow
Reactive Species	Human
Human Gene ID	7145
Mouse Gene ID	21961
Rat Gene ID	301509
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Kidney shows preferential staining in cytoplasm of PCT over DCT. Recommended concentration, 2-4µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

