



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-TERF1 / TRF1 ANTIBODY

SKU: EB09754





SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	telomeric repeat binding protein 1 telomeric repeat binding factor 1 t-TRF1 hTRF1-AS TRF1 TRF TRBF1 PIN2 FLJ41416 telomeric repeat binding factor (NIMA-interacting) 1 TERF1
Accession ID	NP_059523.2
Blocking Peptide	EBP09754
Immunogen	Peptide with sequence C-NLDTRKSVSDKQS, from the internal region of the protein sequence according to NP_059523.2.
Product Comments	This antibody is expected to recognize isoform 1 (NP_059523.2).
Peptide Sequence	C-NLDTRKSVSDKQS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	7013
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 65kDa band observed in Human Testis lysates (calculated MW of 50.2kDa according to NP_059523.2). The observed aberrant MW may be explained by the acetylated form of the protein. Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

