

GOAT ANTI-TERT ANTIBODY

SKU: EB05141



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	telomerase catalytic subunit hEST2 TCS1 EST2 TRT TP2 telomerase reverse transcriptase TERT
Names	
Accession ID	NP_937983.2; NP_937986.1
Immunogen	Peptide with sequence CNPALPSDFKTILD, from the C Terminus of the protein sequence according to NP_937983.2; NP_937986.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_937983.2 ; NP_937986.1).
Peptide Sequence	CNPALPSDFKTILD
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	7015
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Western blot: Approx 130kDa and 75kDa bands seen in HeLa cell lysate. Recommended for use at 1-3 µg/ml.