

Telephone: (650) 697-3600



GOAT ANTI-TFAP2A (ZEBRAFISH) ANTIBODY

SKU: EB07950

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

mont blanc|mob|low|lockjaw|ap2a|ap2|wu:fc31a07|wu:fb83f04|fc31a07|fb83f04|cb93|transcription factor AP-2 **Alias**

alpha|tfap2a **Names**

Accession

NP 789829.1

Blocking

EBP07950 **Peptide**

Peptide with sequence C-NPNSHSETGSKAGDK, from the C Terminus of the protein sequence according to Immunogen

NP 789829.1.

Peptide

C-NPNSHSETGSKAGDK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Zebrafish

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Western **Blot**

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.

Application

Type

Pep-ELISA