

GOAT ANTI-TFAP2A (ZEBRAFISH) ANTIBODY

SKU: EB07950



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	mont blanc mob low lockjaw ap2a ap2 wu:fc31a07 wu:fb83f04 fc31a07 fb83f04 cb93 transcription factor AP-2 alpha tfap2a
Accession ID	NP_789829.1
Blocking Peptide	EBP07950
Immunogen	Peptide with sequence C-NPNSHSETGSKAGDK, from the C Terminus of the protein sequence according to NP_789829.1.
Peptide Sequence	C-NPNSHSETGSKAGDK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Zebrafish
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)