

GOAT ANTI-TFAP2D (AA70-82) ANTIBODY

SKU: EB08868



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	transcription factor AP-2 beta-like 1 transcription factor AP-2 beta (activating enhancer-binding protein 2 beta)-like 1 transcription factor AP-2 beta (activating enhancer binding protein 2 beta)-like 1 activating enhancer binding protein 2 beta-like 1 OTTHUMP00000016593 AP-2 like TFAP2BL1 transcription factor AP-2 delta (activating enhancer binding protein 2 delta) TFAP2D
Accession ID	NP_758438.2
Blocking Peptide	EBP08868
Immunogen	Peptide with sequence C-HHQSFHVEFQHS, from the internal region of the protein sequence according to NP_758438.2.
Peptide Sequence	C-HHQSFHVEFQHS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	83741
Mouse Gene ID	226896
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)