



Telephone: (650) 697-3600

GOAT ANTI-THAP11 / RONIN ANTIBODY

SKU: EB09656

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

THAP domain-containing protein 11|RONIN|HRIHFB2206|CTG-B45d|CTG-B43a|THAP domain containing 11|THAP11 Alias

Names

Accession NP 065190.2 ID

Blocking

EBP09656 **Peptide**

Peptide with sequence C-GSDHSYSLSSGTTEE, from the internal region of the protein sequence according to Immunogen

NP 065190.2.

Peptide

C-GSDHSYSLSSGTTEE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Cow **Species**

Human

57215 Gene ID

Mouse

59016

Gene ID Rat Gene ID 307806

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:1000.

Limit

Application Pep-ELISA

Type