

GOAT ANTI-THAP11 / RONIN ANTIBODY

SKU: EB09656



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	THAP domain-containing protein 11 RONIN HRIHFB2206 CTG-B45d CTG-B43a THAP domain containing 11 THAP11
Accession ID	NP_065190.2
Blocking Peptide	EBP09656
Immunogen	Peptide with sequence C-GSDHSYSLSSGTTEE, from the internal region of the protein sequence according to NP_065190.2.
Peptide Sequence	C-GSDHSYSLSSGTTEE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	57215
Mouse Gene ID	59016
Rat Gene ID	307806
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)