

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600


GOAT ANTI-THRB (ZEBRAFISH) ANTIBODY

SKU: EB07966

Adult hearts samples

 $50 \, \mathrm{kd}$

37 kd



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

wu:fj24f03|trb|fj24f03|TR[b]1|TR-beta|NR1A2|thyroid hormone receptor beta|thrb

Names

Accession

NP_571415.1

Blocking

EBP07966

Peptide

Immunogen

Peptide with sequence C-DIGSAPIVNAPEGNKVD, from the internal region of the protein sequence according to

NP 571415.1.

Product

Comments

This antibody is not expected to cross-react to Thra.

Peptide Sequence

C-DIGSAPIVNAPEGNKVD

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted Species

Zebrafish

Reactive

Zebrafish

Species Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Approx 50kDa band observed in Zebrafish adult Heart lysates (calculated MW of 44.3kDa according to

Blot

NP_571415.1). Recommended concentration: 2-5 μ g/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES

Adult hearts samples

50 kd 37 kd

