


GOAT ANTI-THRB (ZEBRAFISH) ANTIBODY

SKU: EB07966

Adult hearts
samples

50 kd

37 kd



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias wu:fj24f03|trb|fj24f03|TR[b]1|TR-beta|NR1A2|thyroid hormone receptor beta|thrb

Names

Accession ID NP_571415.1

Blocking Peptide EBP07966

Immunogen Peptide with sequence C-DIGSAPIVNAPEGNKVD, from the internal region of the protein sequence according to NP_571415.1.

Product Comments This antibody is not expected to cross-react to Thra.

Peptide Sequence C-DIGSAPIVNAPEGNKVD

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Zebrafish
Reactive Species	Zebrafish
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 50kDa band observed in Zebrafish adult Heart lysates (calculated MW of 44.3kDa according to NP_571415.1). Recommended concentration: 2-5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

Adult hearts
samples

