




www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-THROMBOSPONDIN (INTERNAL) ANTIBODY

SKU: EB08064





SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	thrombospondin-1p180 TSP1 TSP THBS thrombospondin 1
Accession ID	NP_003237.2
Blocking Peptide	EBP08064
Immunogen	Peptide with sequence EGEARETKACKKD, from the internal region of the protein sequence according to NP_003237.2.
Peptide Sequence	EGEARETKACKKD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	7057
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	In transfected HEK293 transiently expressing Human THBS1 a band of approx. 200kDa is observed. This band is not observed in the non-transfected HEK293. The calculated molecular size is 129kDa according to NP_003237.2. Recommended concentration, 1-3µg/ml.
Application Type	Pep-ELISA, WB-trf

DOCUMENTS

- [Data Sheet](#)



GALLERY IMAGES

