

# GOAT ANTI-THYROID PEROXIDASE, BIOTINYLATED ANTIBODY

**SKU:** EB11068-B



250kDa  
150kDa  
100kDa  
75kDa  
50kDa  
37kDa  
25kDa  
20kDa  
15kDa

## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	thyroperoxidase thyroid microsomal antigen TPX TDH2A MSA thyroid peroxidase TPO
<b>Accession ID</b>	NP_000538.3; NP_783650.1; NP_783652.1; NP_783653.1
<b>Blocking Peptide</b>	EBP11068-B
<b>Immunogen</b>	Peptide with sequence ELRCGKHQAVGTSPQ., from the C Terminus of the protein sequence according to NP_000538.3; NP_783650.1; NP_783652.1; NP_783653.1.
<b>Product Comments</b>	This antibody is expected to recognize all reported isoforms (NP_000538.3; NP_783650.1; NP_783652.1; NP_783653.1). Reported variants represent identical protein: NP_000538.3, NP_001193673.1 Reported variants represent identical protein: NP_783650.1, NP_001193674.1
<b>Peptide Sequence</b>	ELRCGKHQAVGTSPQ.
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Reactive Species</b>	Human
<b>Human Gene ID</b>	7173
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:32000.
<b>Western Blot</b>	Approx 110kDa band observed in Human Thyroid Gland lysates (calculated MW of 103kDa according to NP_000538.3). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 1-3µg/ml.
<b>Application Type</b>	Pep-ELISA, WB

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

