

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TIAL1 & TIA ANTIBODY

SKU: EB09793



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

Telephone: (650) 697-3600



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

TIAL1|TIAR|TIA-1-related nucleolysin|TIA-1 related protein|TIA1 cytotoxic granule-associated RNA-binding

Alias p|TIA1 cytotoxic granule-associated RNA binding protein-like 1|T-cluster binding protein|TCBP|OTTHUMP00000020605|OTTHUMP00000020602|MGC33401 Names

Accession

NP_003243.1; NP_001029097.1

Blocking

ID

EBP09793 **Peptide**

Peptide with sequence C-SPEITTEDIKSAFAP, from the internal region of the protein sequence according to **Immunogen**

NP_003243.1; NP_001029097.1.

Product Comments

This antibody is expected to recognize both isoforms (NP 003243.1; NP 001029097.1).

Peptide

Sequence

C-SPEITTEDIKSAFAP

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Dog, Cow

Reactive

Species

Human Gene ID

7073, 7072

Human

Mouse

21843 Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

In paraffin embedded Human Skin shows staining of select keratinocytes. Recommended concentration, **IHC Results**

3-5µg/ml.

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Western Approx 45kDa band observed in lysates of cell line K562 (calculated MW of 43.4kDa according to

Blot NP_001029097.1). Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB, IHC **Type**



Telephone: (650) 697-3600



GALLERY IMAGES

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

