

GOAT ANTI-TICAM-1 / TRIF ANTIBODY

SKU: EB12431



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names proline-rich, vinculin and TI|TIR domain-containing adapter protein inducing IFN-beta|TIR domain-containing adapter molecule 1|TIR domain containing adaptor inducing interferon-beta|TRIF|TICAM-1|PRVTIRB|toll-like receptor adaptor molecule 1|TICAM1

Accession ID NP_891549.1

Blocking Peptide EBP12431

Immunogen Peptide with sequence PDGATFCEDFQVP, from the internal region of the protein sequence according to NP_891549.1.

Peptide Sequence PDGATFCEDFQVP

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog, Pig, Cow

Reactive Species Human

Human Gene ID 148022

Mouse Gene ID 106759

Rat Gene ID 363328

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:128000.

Western Blot Approx 70kDa band observed in lysates of cell line Daudi (calculated MW of 76.4kDa according to NP_891549.1). Recommended concentration: 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

