


GOAT ANTI-TICAM1 ANTIBODY

SKU: EB09579



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TIR domain containing adaptor inducing interferon-beta TRIF TICAM-1 PRVTIRB MGC35334 toll-like receptor adaptor molecule 1 TICAM1
Usage Summary	<p>Flow Cytometry: Flow cytometric analysis of K562 cells. Recommended concentration: 10ug/ml.</p>
Accession ID	NP_891549.1
Blocking Peptide	EBP09579
Immunogen	Peptide with sequence C-HARADEHIALRVREK, from the internal region of the protein sequence according to NP_891549.1.

Peptide Sequence	C-HARADEHIALRVREK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	148022
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png
IHC Results	Paraffin embedded Human Small Intestine. Recommended concentration: 2.5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, FC, IHC

GALLERY IMAGES

