



Telephone: (650) 697-3600



GOAT ANTI-TIM3 ANTIBODY

SKU: EB05339

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

hepatitis A virus cellular receptor 2|kidney injury molecule-3|FLJ14428|KIM-3|T cell immunoglobulin **Alias**

mucin-3|TIM3 **Names**

Accession

NP_116171.3

Peptide with sequence QQPSQPLGCRFAMP, from the C Terminus of the protein sequence according to **Immunogen**

NP 116171.3.

Peptide

QQPSQPLGCRFAMP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Species

Human

Human

84868

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

Western

Western blot: Approx 30-35kDa band seen in Human Brain lysate. Recommended for use at 1-3µg/ml. **Blot**