

GOAT ANTI-TIM3 ANTIBODY

SKU: EB05339



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	hepatitis A virus cellular receptor 2 kidney injury molecule-3 FLJ14428 KIM-3 T cell immunoglobulin mucin-3 TIM3
Accession ID	NP_116171.3
Immunogen	Peptide with sequence QQPSQPLGCRFAMP, from the C Terminus of the protein sequence according to NP_116171.3.
Peptide Sequence	QQPSQPLGCRFAMP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	84868
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 30-35kDa band seen in Human Brain lysate. Recommended for use at 1-3µg/ml.