

GOAT ANTI-TIM3 (MOUSE) ANTIBODY

SKU: EB05340



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	T-cell immunoglobulin and mucin domain containing 3 hepatitis A virus cellular receptor 2 Timd3 Tim3 Havcr2
Accession ID	NP_599011
Immunogen	Peptide with sequence NSNEYCYVNSQQPS, from the C Terminus of the protein sequence according to NP_599011.
Peptide Sequence	NSNEYCYVNSQQPS
Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat
Mouse Gene ID	171285
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Jurkat and MOLT-4 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?