


GOAT ANTI-TINF2 ANTIBODY

SKU: EB06920



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names TRF1-interacting nuclear protein 2| TERF1-interacting nuclear factor 2| (TRF1)-interacting nuclear factor 2 variant 1| TIN2| DKCA3| TERF1 (TRF1)-interacting nuclear factor 2|TINF2

Accession ID NP_036593.1

Blocking Peptide EBP06920

Immunogen Peptide with sequence C-KENPVDLPATEQKE, from the C Terminus of the protein sequence according to NP_036593.1.

Peptide Sequence C-KENPVDLPATEQKE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	26277
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Heart and Kidney. Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

