


GOAT ANTI-TLR4 / CD284 ANTIBODY

SKU: EB09441



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	homolog of Drosophila toll OTTHUMP00000022807 hToll TOLL CD284 ARMD10 toll-like receptor 4 TLR4
Usage Summary	Immunofluorescence: Strong expression of the protein seen in the cytoplasm of NIH3T3 cells and membrane and cytoplasm of THP-1 cells. Recommended concentration: 10µg/ml Flow Cytometry: Flow cytometric analysis of U937 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_612564.1
Blocking Peptide	EBP09441
Immunogen	Peptide with sequence C-YSSQDEDWVRNE, from the Internal region of the protein sequence according to NP_612564.1.
Peptide Sequence	C-YSSQDEDWVRNE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human, Mouse
Human Gene ID	7099
Mouse Gene ID	21898
Rat Gene ID	29260
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png
IHC Results	Paraffin embedded Human Spleen. Recommended concentration: 5-6µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IF, IHC, FC

GALLERY IMAGES

